


Omega Optical Filter Set XF06



FILTER SET COMPONENTS

Exciter: XF1005 (365WB50)

Dichroic: XF2001 (400DCLP)

Emitter: XF3002 (450DF65)

OTHER INFORMATION

Bandpass blue set with a 365nm wide-band exciter. Higher signal-to-noise than the XF05-2. Highly recommended for multicolour work using a mercury (Hg) light source. This set will also work for a xenon (Xe) light source, but in this case, the set will have a much dimmer signal than our XF03 filter set, which is more optimised for xenon arc lamps. For use with Hg arc lamp.