


## Omega Optical Filter Set XF104-3



### FILTER SET COMPONENTS

Exciter:	XF1068 (500AF25)
Dichroic:	XF2030 (525DRLP)
Emitter:	XF3106 (545AF35)

### OTHER INFORMATION

**Notes:** This set XF104-2 option incorporates ALPHA technology. The set has a better signal to noise ratio than our previous XF104 set. This set is excellent for work with YFP (yellow fluorescent protein) and other dyes listed above. This set will also work well for multilabeling with DsRed, the new red fluorescent protein from Clontech Laboratories. In addition, this set provides excellent spectral discrimination from CFP (cyan fluorescent protein). If you prefer a longpass set to maximise signal, see XF105-2.