


## Omega Optical Filter Set XF108-2



### FILTER SET COMPONENTS

Exciter:	XF1074 (525AF45)
Dichroic:	XF2017 (560DRLP)
Emitter:	XF3083 (595AF60)

### OTHER INFORMATION

Highest quality orange bandpass set for most rhodamine fluorophores. This set XF108-2 option incorporates ALPHA technology. The set has a better signal to noise ratio than our previous XF108 set. This set is a bandpass alternative to our XF101-2 offering. This set is recommended for multilabeling situations in which one might want to use Cy5 or other dyes to the red of Cy3. Emission colour is orange-yellow.