


Omega Optical Filter Set XF111-2



FILTER SET COMPONENTS

Exciter: XF1077 (540AF30)

Dichroic: XF2015 (570DRLP)

Emitter: XF3089 (575ALP)

OTHER INFORMATION

This set is designed for imaging DsRed (RFP) with a black and white detector, such as a CCD, film, video camera, or PMT. It will also adequately separate DsRed signal from other FPs. The emission colour of this set is orange-yellow. This set incorporates our new ALPHA technology, resulting in a superior signal to noise ratio.