


Omega Optical Filter Set XF114-2



FILTER SET COMPONENTS

Exciter:	XF1071 (440AF21)
Dichroic:	XF2034 (455DRLP)
Emitter:	XF3075 (480AF30)

OTHER INFORMATION

Highest quality CFP bandpass set. Designed for superb YFP discrimination. This set XF114-2 option incorporates ALPHA technology. The set has a better signal to noise ratio than our previous XF114 set. This set features a pale blue signal with a very black background. It provides excellent spectral separation from YFP. See XF130 for a CFP longpass filter set.