


Omega Optical Filter Set XF124



<u>FILTER SET COMPONENTS</u>	
Exciter:	XF1006 (400DF15)
Dichroic:	XF2059 (415DRLP)
Emitter:	XF3002 (450DF65)
Emitter2:	XF3007 (535DF35)

OTHER INFORMATION

Dual emission set designed to discriminate blue and green of GeneBLAzer, a beta lactamase assay developed by Dr. Roger Tsien of UCSD. See XF106 for a set that will allow simultaneous visualisation of the blue and green emissions.