


## Omega Optical Filter Set XF154-1



<b>FILTER SET COMPONENTS</b>	
Exciter:	XF1201 (436AF8)
Exciter:	XF1042 (485DF15)
Exciter:	XF 1044 (575DF25)
Dichroic:	XF2090 (445-510-600TBDR)
Emitter:	XF3118 (465-535-640TBEM)

### **OTHER INFORMATION**

A series of singleband exciters in a filter wheel or slider with a multi-band dichroic & emission filter housed in the microscope. This set is used to collect images for CFP-YFP-DsRed2 with a PMT or CCD camera. Letter designations are custom filter specific for the filter set being shipped & may not necessarily match the spectral characterisations exactly as displayed. Please note these exciters can not be ordered individually.