


Omega Optical Filter Set XF300



FILTER SET COMPONENTS

Exciter1:	XF1301 (415WB100)
Exciter2:	XF1009 (425DF45)
Dichroic:	XF2007 (475DCLP)
Emitter:	XF3301 (525WB20)
Emitter:	XF3002 (565WB20)
Emitter:	XF3003 (585WB20)
Emitter:	XF3004 (605WB20)
Emitter:	XF3005 (665WB20)
Emitter:	XF3113 (710AF40)
Emitter:	XF3307 (800WB80)
Emitter:	XF3308 (840WB80)

OTHER INFORMATION

Set can be used for Qdot™ 525, 565, 585, 605, 655, 705, & 800 Conjugates.

