


Omega Optical Filter Set XF53-1



FILTER SET COMPONENTS

Exciter:	XF1042 (485DF15)
Exciter2:	XF1044 (575DF25)
Dichroic:	XF2044 (490-575DBDR)
Emitter:	XF3057 (528-633DBEM)

OTHER INFORMATION

A series of singleband exciters in a filter wheel or slider with a multi-band dichroic & emission filter housed in the microscope. This set is used to collect images for FITC-Texas Red with a PMT or a CCD camera. Letter designations are custom filter specific for the filter set being shipped & may not necessarily match the spectral characteristics exactly as displayed. Please note these exciters can not be ordered individually..