


Omega Optical Filter Set XF63-1



FILTER SET COMPONENTS

Exciter: XF1006 (400DF15)

Exciter2: XF1042 (485DF15)

Exciter3: XF1044 (575DF25)

Dichroic: XF2048 (400-477-575TBDR)

Emitter: XF3061 (445-525-650TBEM)

OTHER INFORMATION

Pinkel Set for DAPI/FITC/Texas Red. A series of singleband excitors in a filter wheel or slider with a multi-band dichroic & emission filter housed in the microscope. This set is used to collect images for DAPI-FITC-Texas Red with a PMT or a CCD camera.

Letter designations are custom filter specific for the filter set being shipped & may not necessarily match the spectral characteristics exactly as displayed.
Please note these excitors can not be ordered individually.