


## Omega Optical Filter Set XF67-1



### FILTER SET COMPONENTS

Exciter:	XF1006 (400DF15)
Exciter2:	XF1042 (485DF15)
Exciter3:	XF1044 (575DF25)
Dichroic:	XF2045 (400-485-580TBDR)
Emitter:	XF3058 (457-528-633TBEM)

### OTHER INFORMATION

Pinkel Set for DAPI/FITC/Texas Red. A series of singleband exciters in a filter wheel or slider with a multi-band dichroic & emission filter housed in the microscope. This set is used to collect images for DAPI-FITC-Texas Red, DAPI-Alexa Fluor 488/546 or DAPI-Cy2-Cy3 with a PMT or a CCD. Letter designations are custom filter specific for the filter set being shipped & may not necessarily match the spectral characteristics exactly as displayed. Please note these exciters can not be ordered individually.