


# Omega Optical Filter Set XF92



| <b><u>FILTER SET COMPONENTS</u></b> |                      |
|-------------------------------------|----------------------|
| <b>Exciter:</b>                     | XF1062 (550-640DBEX) |
| <b>Dichroic:</b>                    | XF2053 (555-640DBDR) |
| <b>Emitter:</b>                     | XF3066 (595-700DBEM) |

## **OTHER INFORMATION**

This set is designed to allow simultaneous imaging of Cy3 and Cy5. Please be advised that this set works best with a CCD detector since the emission of Cy5 is too far to the red for the eye to see.