

# PHOTO RESEARCH®

## Camera Kit



Camera Kit Shown on PR-705

The Photo Research Camera Kit provides an easy method for remotely monitoring the field of view of your Photo Research photometer or spectroradiometer.

The adapter component of the system easily marries the "lipstick" CCD camera to the instrument eyepiece.

Simply hook up the camera connector to the system monitor via the standard BNC connector and begin monitoring the scene. Alternately, hook up the camera output to the video board in your PC (video board not included) to capture the measurement scene on your desktop or laptop computer. This can be a powerful tool for correlating measurements between different locations by illustrating precisely what was measured.

### The Camera Kit Includes:

- Eyepiece Adapter
- CCD Camera (B/W or Color)
- Camera Power Supply
- Monitor (B/W or Color)
- All Connecting Cables
- Instructions

### The Camera Kit is Available For:

- PR-650
- PR-703/713
- PR-704/714
- PR-705/715
- PR-880
- PR-1500/1530AR/1600
- PR-1980A/1980B/1980C



Camera Kit with Live Image  
on System Monitor  
(camera power supply not shown)



**PHOTO RESEARCH, INC.**

9731 Topanga Canyon Place Chatsworth, CA 91311-4125 USA  
Phone: 818/341-5151 • Fax: 818/341-7070  
<http://www.photoresearch.com>  
e-mail: [sales@photoresearch.com](mailto:sales@photoresearch.com)